









LOCKING ARM MOVES
ASIDE TO ALLOW THE
SCREW TO PASS, THEN
"CLICKS" BACK OVER TO
CAPTURE AND LOCK THE
SCREW.

SIZES

MEASURED END TO END

LEVELS	SIZE RANGE	INCREMENTS
1 LEVEL	22-34mm	2mm
2 LEVEL	36-54mm	3mm
3 LEVEL	53-77mm	3mm
4 LEVEL	68-92mm	4mm

SCREWS

SELF TAPPING & SELF DRILLING LENGTHS AVAIABLE:

12mm, 14mm, 16mm (18mm UPON REQUEST)



PLATES		
Item #	Description	
PSN60000100	22mm, Single Level Plate	
PSN60000101	24mm, Single Level Plate	
PSN60000102	26mm, Single Level Plate	
PSN60000103	28mm, Single Level Plate	
PSN60000104	30mm, Single Level Plate	
PSN60000105	32mm, Single Level Plate	
PSN60000106	34mm, Single Level Plate	
PSN60000107	36mm, Two Level Plate	
PSN60000108	39mm, Two Level Plate	
PSN60000109	42mm, Two Level Plate	
PSN60000110	45mm, Two Level Plate	
PSN60000111	48mm, Two Level Plate	
PSN60000112	51mm, Two Level Plate	
PSN60000113	54mm, Two Level Plate	
PSN60000114	53mm, Three Level Plate	
PSN60000115	56mm, Three Level Plate	
PSN60000116	59mm, Three Level Plate	
PSN60000117	62mm, Three Level Plate	
PSN60000118	65mm, Three Level Plate	
PSN60000129	68mm, Three Level Plate	
PSN60000119	71mm, Three Level Plate	
PSN60000120	74mm, Three Level Plate	
PSN60000121	77mm, Three Level Plate	
PSN60000122	68mm, Four Level Plate	
PSN60000123	72mm, Four Level Plate	
PSN60000124	76mm, Four Level Plate	
PSN60000125	80mm, Four Level Plate	
PSN60000126	84mm, Four Level Plate	
PSN60000127	88mm, Four Level Plate	
PSN60000128	92mm, Four Level Plate	

SCREWS		
ltem #	Description	
PSN60000135	4 X 12mm Var Screw	
PSN60000136	4 X 14mm Var Screw	
PSN60000137	4 x 16mm Var Screw	
PSN60000140	4.25 x 12mm Var Screw	
PSN60000141	4.25 x 14mm Var Screw	
PSN60000142	4.25 x 16mm Var Screw	
PSN60000145	4 X 12mm Fixed Screw	
PSN60000146	4 X 14mm Fixed Screw	
PSN60000147	4 x 16mm Fixed Screw	
PSN60000150	4.25 x 12mm Fixed Screw	
PSN60000151	4.25 x 14mm Fixed Screw	
PSN60000152	4.25 x 16mm Fixed Screw	
PSN60001031	4 x 18mm Var Screw	
PSN60001032	4.25 x 18mm Var Screw	
PSN60001033	4 X 18mm Fixed Screw	
PSN60001034	4.25 x 18mm Fixed Screw	



Presidio Surgical, Inc. 1451 Danville Blvd. Suite 204 Alamo, CA 94507

PHONE: (916)941-1967 FAX: (925)217-4821

www.presidiosurgicalspine.com

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Products Patented and Patents Pending

All products are not currently available in all markets